

App Guide  
**Standing Balance**  
**Trunk Front and Side Angle**  
SageMotion  
Wearable Biofeedback System





## **Table of Contents**

Components	<i>(page 1)</i>
Wirelessly Connect to Computer or Cellphone	<i>(page 2)</i>
Standing Balance Trunk Front and Side Angle	<i>(page 3)</i>

# Components



Hub



Nodes (8x)



Battery



Node Straps: *Medium (8x), Short (4x), Long (2x)*



Cable A (10x)

-Connect Hub to Battery  
-Charge Nodes & Battery



Cable B (*optional use*)

-Connect Hub to Computer




Node Charging Station

# Wirelessly Connect to Computer or Cellphone

## 1) Connect Cable A to Battery and to Hub




## 2) On Computer/Cellphone, Connect to Wi-Fi: "SageMotion"



*Note 1: Need to wait for up to 1 minute for "SageMotion" to appear in Wi-Fi list. If it doesn't appear, try turning the Wi-Fi off and then on again on the computer/cellphone.*

*Note 2: Hub is connected after clicking "Connect" even if in Windows it shows "Connecting" or "No internet, open".*

## 3) On Computer/Cellphone, in Chrome Address Bar, Go To <http://192.168.12.1>



**[Note] If Computer Doesn't Have Wi-Fi:** plug in Cable B to the Hub and to the ethernet port of your computer, then in chrome address bar, go to **<http://192.168.137.1>**

# Standing Balance Trunk Front and Side Angle App

*The purpose of the Standing Balance Trunk Front and Side Angle App is to record, analyze, and provide feedback for Trunk Side Angle, Trunk Front Angle while subjects perform Standing Balance activities.*

## 1) Turn on 5 Nodes




*Slide switch toward middle to turn node on*



*Green light will blink after the node is on and running*

## 2) Select “Standing Balance Trunk Front and Side Angle” App



## 3) Click “Search”

**Node List**



# Standing Balance Trunk Front and Side Angle App (cont.)

4) Configure 1 Sensor Nodes and 4 Feedback Nodes as Shown Below:

Node List

Search

Connect

Type	Position	MAC	
<div>sensor</div>	<div>trunk</div>	88:6B:0F:E1:D8:A2	<div></div>
<div>feedback</div>	<div>feedback_right</div>	88:6B:0F:E1:D8:9E	<div></div>
<div>feedback</div>	<div>feedback_left</div>	88:6B:0F:E1:D8:A6	<div></div>
<div>feedback</div>	<div>feedback_front</div>	88:6B:0F:E1:D8:9F	<div></div>
<div>feedback</div>	<div>feedback_back</div>	88:6B:0F:E1:D8:96	<div></div>

5) Click “Connect”

Node List

Search

Connect

6) “Ready to collect data” Will Appear after Node Connection is Complete

Standing Balance 1

Start



Stop

✓ Ready to collect data

# Standing Balance Trunk Front and Side Angle App (cont.)





7) Thread Straps through Nodes and Attach at Locations Shown Below:

*Trunk sensor orientations on/off switch must point upwards*





















*trunk node placed on the back of the torso*

### How to Thread Straps



[Note] feedback\_left, feedback\_right, feedback\_front, feedback\_back nodes can be placed at any location

8) Click “Blink” for each Node to Confirm Correct Locations (red LED for given node blinks 3 times on click)

Type	Position	MAC			
sensor	trunk	88:6B:0F:E1:D8:A2			
feedback	feedback_right	88:6B:0F:E1:D8:9E			
feedback	feedback_left	88:6B:0F:E1:D8:A6			
feedback	feedback_front	88:6B:0F:E1:D8:9F			
feedback	feedback_back	88:6B:0F:E1:D8:96			



# Standing Balance Trunk Front and Side Angle App (cont.)

9) In App Configuration, Enter Settings (Example Below)

## App Configuration

Trial Name	<input type="text" value="trial_1"/>
<b>Feedback Settings</b>	
Feedback On	<input type="text" value="true"/>
Max Lean Right Angle	<input type="text" value="10"/>
Max Lean Left Angle	<input type="text" value="10"/>
Max Lean Front Angle	<input type="text" value="10"/>
Max Lean Back Angle	<input type="text" value="10"/>
<b>Save Options</b>	
Save Mode	<input type="text" value="xlsx"/>

# Standing Balance Trunk Front and Side Angle App (cont.)

10) Click “Start” to Start Running the App



11) After the Trial is Finished, Click “Stop”



12) After Clicking “Stop”, a File from that Trial will Appear under Download Data. Click the File (e.g. trial\_1) to Download it to the Computer or Phone.

Data List									
<input type="checkbox"/> Name	Date▲	Duration	App	Type	Size	Rename	Delete		
<input type="checkbox"/> Trial 1	2024-07-12-13-29-39	0:00:06	Standing Balance Trunk Front and Side Angle	.xlsx	619.8 kB				

# Standing Balance Trunk Front and Side Angle App (*cont.*)

## Description of Data in Downloaded File

**time** (sec): time since trial start

**TSA** (deg): trunk side angle (medial-lateral), positive is to the right

**TFA** (deg): trunk forward angle (anterior-posterior), positive is forward

**max\_lean\_right** (deg): The feedback threshold of max right lean

**max\_lean\_left** (deg): The feedback threshold of max left lean

**max\_lean\_front** (deg): The feedback threshold of max front lean

**max\_lean\_back** (deg): The feedback threshold of max back lean

**SensorIndex\_1**: index of raw sensor data

**AccelX/Y/Z\_1** (m/s<sup>2</sup>): raw acceleration data

**GyroX/Y/Z\_1** (deg/s): raw gyroscope data

**MagX/Y/Z\_1** (μT): raw magnetometer data

**Quat1/2/3/4\_1**: quaternion data

**Sampletime\_1**: timestamp of the sensor

**Package\_1**: package number of the sensor